

## Histone H3 Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP10433
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Purified recombinant Full length Histone H3 (Human)
<b>Mol wt</b>	N/A
<b>Species reactivity</b>	Human, Mouse, Rat, Yeast
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, IP
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Histone H3
<b>Synonyms</b>	Histone H3; Histone H3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

#### Recommended Dilution


WB: 1: 5000 - 1: 20000

IP: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

### Images



Western Blot analysis of various cells using Histone H3 rabbit pAb. Secondary antibody was diluted at 1:20000

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)