

Histone H3 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10433
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Purified recombinant Full length Histone H3 (Human)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat, Yeast
Clonality	Polyclonal
Recommended application	WB, IP
Concentration	1 mg/ml
Full name	Histone H3
Synonyms	Histone H3; Histone H3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution


WB: 1: 5000 - 1: 20000

IP: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of various cells using Histone H3 rabbit pAb. Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com