

IL1RA Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10584
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human IL1RA
Mol wt	19470
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	IL1RA
Synonyms	IL1RA; Interleukin-1 receptor antagonist protein; IL-1RN; IL-1ra; IRAP; ICIL-1RA; IL1 inhibitor; Anakinra

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Genetic variation in IL1RN is associated with susceptibility to diabetic nephropathy

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com