

# IQEC1 Rabbit Polyclonal Antibody

## Description

Product type	Primary Antibody
Code	BT-AP10660
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human IQEC1
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	IQEC1
Synonyms	IQEC1; IQ motif and SEC7 domain-containing protein 1; ADP-ribosylation factors guanine nucleotide-exchange protein 100; ADP-ribosylation factors guanine nucleotide-exchange protein 2; Brefeldin-resistant Arf-GEF 2 protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.


## Background

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

## Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4°C overnight

## Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)