

Jak2(Phospho-Tyr1008) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10725
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human Jak2 (Tyr1008)
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Jak2
Synonyms	Jak2 ;Tyr1008 ; Tyrosine-protein kinase JAK2; EC 2.7.10.2; Janus kinase 2; JAK-2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com