

K1161 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10740
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human K1161
Mol wt	78540
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Uncharacterized family 31 glucosidase KIAA1161
Synonyms	Uncharacterized family 31 glucosidase KIAA1161 ;EC 3.2.1.-

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year