

KIF1B(Phospho-Ser1487) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10829
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human KIF1B (Ser1487)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	KIF1B
Synonyms	KIF1B ;Ser1487; Kinesin-like protein KIF1B; Klp

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a motor protein that transports mitochondria and synaptic vesicle precursors. Mutations in this gene cause Charcot-Marie-Tooth disease| type 2A1.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com