

LRP1B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP10972
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 200-280
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF
Concentration	1 mg/ml
Full name	Low-density lipoprotein receptor-related protein 1B
Synonyms	Low-density lipoprotein receptor-related protein 1B ;LRP-1B;Low-density lipoprotein receptor-related protein-deleted in tumor;LRP-DIT

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a member of the low density lipoprotein (LDL) receptor family. These receptors play a wide variety of roles in normal cell function and development due to their interactions with multiple ligands. Disruption of this gene has been reported in several types of cancer.

Recommended Dilution

IHC-p: 1: 50 - 1: 300

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com