

LIPA4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11049
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human LIPA4
Mol wt	77110
Species reactivity	Human, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	LIPA4
Synonyms	LIPA4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


PPFIA4| or liprin-alpha-4| belongs to the liprin-alpha gene family. See liprin-alpha-1 (LIP1| or PPFIA1; MIM 611054) for background on liprins.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com