

Lyn(Phospho-Tyr507) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11177
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phosho peptide around human Lyn (Tyr507)
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Lyn
Synonyms	Lyn ;Tyr507; Tyrosine-protein kinase Lyn; EC 2.7.10.2; Lck/Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a tyrosine protein kinase| which maybe involved in the regulation of mast cell degranulation| and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year