

MA1A2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11203
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MA1A2
Mol wt	70510
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Mannosyl-oligosaccharide 1 2-alpha-mannosidase IB
Synonyms	Mannosyl-oligosaccharide 1; 2-alpha-mannosidase IB ;EC 3.2.1.113;Mannosidase alpha class 1A member 2;Processing alpha-1; 2-mannosidase IB;Alpha-1; 2-mannosidase IB

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com