

MERIT40(Phospho-Ser29) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11290
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human MERIT40 (Ser29)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	MERIT40
Synonyms	MERIT40 ;Ser29; BRISC and BRCA1-A complex member 1; Mediator of RAP80 interactions and targeting subunit of 40 kDa; New component of the BRCA1-A complex

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com