

MET2B Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11299
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MET2B
Mol wt	41580
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	MET2B
Synonyms	MET2B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene is a member of a family of methyltransferases that share homology with| but are distinct from| the UbiE family of methyltransferases. Alternatively spliced variants which encode different protein isoforms have been described; however| not all variants have been fully characterized.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from A549 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com