

MRP4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11370
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MRP4
Mol wt	145750
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Multidrug resistance-associated protein 4
Synonyms	Multidrug resistance-associated protein 4 ;ATP-binding cassette sub-family C member 4;MRP/cMOAT-related ABC transporter;Multi-specific organic anion transporter B;MOAT-B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaying, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com