

MRT4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11371
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MRT4
Mol wt	26290
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	MRT4
Synonyms	MRT4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene encodes a protein sharing a low level of sequence similarity with ribosomal protein P0. While the precise function of the encoded protein is currently unknown it appears to be involved in mRNA turnover and ribosome assembly.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com