

MUC6 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11424
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human MUC6
Mol wt	268290
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	MUC6
Synonyms	MUC6; Mucin-6; MUC-6; Gastric mucin-6

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The MUC6 gene encodes gastric mucin| a secreted glycoprotein that plays an essential role in epithelial cytoprotection from acid| proteases| pathogenic microorganisms| and mechanical trauma in the gastrointestinal tract (summary by Toribara et al.| 1993

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com