

MOT11 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11441
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 1-80
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Monocarboxylate transporter 11
Synonyms	Monocarboxylate transporter 11 ;MCT 11;Solute carrier family 16 member 11

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

It is uncertain whether Met-1 or Met-25 is the initiator.,Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates.,Belongs to the major facilitator superfamily. Monocarboxylate porter (TC 2.A.1.13) family.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com