

NACA2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11486
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human NACA2
Mol wt	23650
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Nascent polypeptide-associated complex subunit alpha-2
Synonyms	Nascent polypeptide-associated complex subunit alpha-2 ;Alpha-NAC-like;Hom s 2.01;Nascent polypeptide-associated complex subunit alpha-like;NAC-alpha-like

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com