

NPRL3 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11572
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human NPRL3
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	NPRL3
Synonyms	NPRL3; Nitrogen permease regulator 3-like protein; -14 gene protein; Alpha-globin regulatory element-containing gene protein; Protein CGTHBA

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The function of the encoded protein is not known.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com