

MYBPH Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11600
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Myosin-binding protein H
Synonyms	Myosin-binding protein H ;MyBP-H;H-protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Binds to myosin; probably involved in interaction with thick myofilaments in the A-band.,Belongs to the immunoglobulin superfamily. MyBP family.,Contains 2 fibronectin type-III domains.,Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Skeletal muscle.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com