

# Myc-Tag Polyclonal Antibody, HRP Conjugated

### Description

Product type	Primary Antibody
Code	BT-AP11602
Host	Rabbit
Isotype	IgG
Size	100u1
Immunogen	N/A
Mol wt	N/A
Species reactivity	Species independent
Clonality	Polyclonal
Recommended application	WB
Concentration	1mg/ml
Full name	Myc-Tag
Synonyms	Myc-Tag; Myc

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

# Background

#### **Recommended Dilution**

WB: 1: 5000 ELISA: 1: 20000 Not yet tested in other applications.

## Images

No images.

Storage -20°C for 1 year

> 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com