

# O11H2 Rabbit Polyclonal Antibody

## Description

Product type Primary Antibody

Code BT-AP11617

Host Rabbit

Isotype IgG

Size 100ul, 50ul, 20ul

Immunogen Synthesized peptide derived from human O11H2

Mol wt 35860

Species reactivity Human, Rat, Mouse

**Clonality** Polyclonal

Recommended application WB

Concentration 1 mg/ml

Full name Olfactory receptor 11H2

Synonyms Olfactory receptor 11H2 ;Olfactory receptor OR14-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

## Background

## **Recommended Dilution**

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

#### Images

No images.

#### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com