

NAC2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11687
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 380-460
Mol wt	N/A
Species reactivity	Human, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Sodium/calcium exchanger 2
Synonyms	Sodium/calcium exchanger 2 ;Na ⁺ /Ca ²⁺ -exchange protein 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

enzyme regulation:By ATP.,Rapidly transports Ca²⁺ during excitation-contraction coupling. Ca²⁺ is extruded from the cell during relaxation so as to prevent overloading of intracellular stores.,Belongs to the sodium/potassium/calcium exchanger family. SLC8 subfamily.,Contains 2 Calx-beta domains.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com