

PEO1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11877
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PEO1
Mol wt	75240
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Twinkle protein, mitochondrial
Synonyms	Twinkle protein; mitochondrial ;EC 3.6.4.12;Progressive external ophthalmoplegia 1 protein;T7 gp4-like protein with intramitochondrial nucleoid localization;T7-like mitochondrial DNA helicase

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com