

## PFD4 Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP11890
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PFD4
Mol wt	14740
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	PFD4
Synonyms	PFD4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin| a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides| thereby allowing them to fold correctly. The complex| consisting of two alpha and four beta subunits| forms a double beta barrel assembly with six protruding coiled-coils.

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4°C overnight

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)