

PGCA (Cleaved-Ala17) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP11896
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human PGCA (Cleaved-Ala17)
Mol wt	265650
Species reactivity	Human, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	PGCA
Synonyms	PGCA ;Cleaved-Ala17; Aggrecan core protein; Cartilage-specific proteoglycan core protein; CSPCP; Chondroitin sulfate proteoglycan core protein 1; Chondroitin sulfate proteoglycan 1; Aggrecan core protein 2;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

alternative products:Additional isoforms seem to exist|developmental stage:Expression was detected in chondrocytes throughout the developing skeleton.|disease:Defects in ACAN are the cause of spondyloepiphyseal dysplasia type Kimberley (SEDK)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com