

PPIG(Phospho-Ser376) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12084
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human PPIG (Ser376)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	PPIG
Synonyms	PPIG ;Ser376; Peptidyl-prolyl cis-trans isomerase G; PPIase G; Peptidyl-prolyl isomerase G; EC 5.2.1.8; CASP10; Clk-associating RS-cyclophilin; CARS-Cyp; CARS-cyclophilin; SR-cyclophilin; SR-cyp; SRcyp; Cyclophilin G; Rotamase G

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year