

QCR10 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12220
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human QCR10
Mol wt	6160
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Cytochrome b-c1 complex subunit 10
Synonyms	Cytochrome b-c1 complex subunit 10 ;Complex III subunit 10;Complex III subunit XI;Ubiquinol-cytochrome c reductase complex 6.4 kDa protein

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com