

ROS1(Phospho-Tyr2274) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12336
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human ROS1 (Tyr2274)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	ROS1
Synonyms	ROS1 ;Tyr2274; Proto-oncogene tyrosine-protein kinase ROS; EC 2.7.10.1; Proto-oncogene c-Ros; Proto-oncogene c-Ros-1; Receptor tyrosine kinase c-ros oncogene 1; c-Ros receptor tyrosine kinase

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This proto-oncogene| highly-expressed in a variety of tumor cell lines| belongs to the sevenless subfamily of tyrosine kinase insulin receptor genes. The protein encoded by this gene is a type I integral membrane protein with tyrosine kinase activity. The protein may function as a growth or differentiation factor receptor.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com