

S15A2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12367
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human S15A2
Mol wt	80190
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Solute carrier family 15 member 2
Synonyms	Solute carrier family 15 member 2 ;Kidney H;+/peptide cotransporter;Oligopeptide transporter; kidney isoform;Peptide transporter 2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year