

## S38A6 Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP12412
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human S38A6
<b>Mol wt</b>	50160
<b>Species reactivity</b>	Human, Rat
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Probable sodium-coupled neutral amino acid transporter 6
<b>Synonyms</b>	Probable sodium-coupled neutral amino acid transporter 6 ;N-system amino acid transporter 1;NAT-1;Na <sup>+</sup> -coupled neutral amino acid transporter 6;Solute carrier family 38 member 6

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

#### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year