

## SIA7B Rabbit Polyclonal Antibody

### Description

Product type	Primary Antibody
Code	BT-AP12656
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human SIA7B
Mol wt	41140
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	SIA7B
Synonyms	SIA7B

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background


ST6GALNAC2 belongs to a family of sialyltransferases that add sialic acids to the nonreducing ends of glycoconjugates. At the cell surface| these modifications have roles in cell-cell and cell-substrate interactions| bacterial adhesion| and protein targeting (Samyn-Petit et al. | 2000

### Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

### Images



### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)