

SNX17 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12703
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human SNX17
Mol wt	51700
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	SNX17
Synonyms	SNX17

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain| which is a phosphoinositide binding domain| and are involved in intracellular trafficking. This protein does not contain a coiled coil region| like some family members| but contains a B41 domain. This protein interacts with the cytoplasmic domain of P-selectin| and may function in the intracellular trafficking of P-selectin. Alternative splicing results in multiple transcript variants encoding different isoforms.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year