

## SNX8 Rabbit Polyclonal Antibody

### Description

|                                |   |
|--------------------------------|---|
| <b>Product type</b>            | Primary Antibody                            |
| <b>Code</b>                    | BT-AP12711                                  |
| <b>Host</b>                    | Rabbit                                      |
| <b>Isotype</b>                 | IgG   |
| <b>Size</b>                    | 100ul, 50ul, 20ul                           |
| <b>Immunogen</b>               | Synthesized peptide derived from human SNX8 |
| <b>Mol wt</b>                  | 51150                                       |
| <b>Species reactivity</b>      | Human, Mouse                                |
| <b>Clonality</b>               | Polyclonal                                  |
| <b>Recommended application</b> | WB  |
| <b>Concentration</b>           | 1 mg/ml                                     |
| <b>Full name</b>               | Sorting nexin-8                             |
| <b>Synonyms</b>                | Sorting nexin-8                             |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

#### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

#### Images

No images.

#### Storage

-20°C for 1 year