

SRY (Acetyl Lys136) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12820
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human SRY (Acetyl Lys136)
Mol wt	22440
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	SRY
Synonyms	SRY ;Acetyl Lys136; Sex-determining region Y protein; Testis-determining factor

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in SRY are a cause of gonadal dysgenesis XY female type (GDXY)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com