

PEAR1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP12923
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 910-990
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Platelet endothelial aggregation receptor 1
Synonyms	Platelet endothelial aggregation receptor 1 ;hPEAR1;Multiple epidermal growth factor-like domains protein 12;Multiple EGF-like domains protein 12

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

PEAR1 is a platelet receptor that signals upon the formation of platelet-platelet contacts independent of platelet activation and secondary to platelet aggregation (Nanda et al., 2005)

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com