

TICRR Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13017
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TICRR
Mol wt	210100
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TICRR
Synonyms	TICRR

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Treslin is involved in the initiation of DNA replication (Kumagai et al. | 2010)

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com