

TRAD1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13152
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human TRAD1
Mol wt	64020
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	TRAD1
Synonyms	TRAD1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The innate immune system confers host defense against viral and microbial infection| and TRAFD1 is a negative feedback regulator that controls excessive immune responses (Sanada et al.| 2008

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com