

TRDMT Rabbit Polyclonal Antibody

Description

| | |
|--------------------------------|--|
| Product type | Primary Antibody |
| Code | BT-AP13161 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 100ul, 50ul, 20ul |
| Immunogen | Synthesized peptide derived from human TRDMT |
| Mol wt | 43010 |
| Species reactivity | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Recommended application | WB |
| Concentration | 1 mg/ml |
| Full name | tRNA |
| Synonyms | tRNA ;cytosine;38-C;5-methyltransferase ;EC 2.1.1.204;DNA ;cytosine-5-methyltransferase-like protein 2;Dnmt2;DNA methyltransferase homolog HsaIIP;DNA MTase homolog HsaIIP;M.HsaIIP;PuMe |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year