

TRDMT Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13161
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human TRDMT
Mol wt	43010
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	tRNA
Synonyms	tRNA ;cytosine;38-C;5-methyltransferase ;EC 2.1.1.204;DNA ;cytosine-5-methyltransferase-like protein 2;Dnmt2;DNA methyltransferase homolog HsaIIP;DNA MTase homolog HsaIIP;M.HsaIIP;PuMe

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com