

PLXA4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13191
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Plexin-A4
Synonyms	Plexin-A4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data. Forms a receptor complex with neuropilin-1 and can propagate semaphorin-3A elicited inhibitory signals into cells and neurons. PTM: Phosphorylated upon DNA damage, probably by ATM or ATR. Belongs to the plexin family. Contains 1 Sema domain. Contains 3 PSI domains. Contains 4 IPT/TIG domains.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com