

Tubulin alpha Rabbit Polyclonal Antibody(Caenorhabditis elegans)

Description

Product type	Primary Antibody
Code	BT-AP13265
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Purified recombinant Full length Tubulin alpha (Caenorhabditis elegans)
Mol wt	N/A
Species reactivity	C.elegans
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Tubulin alpha
Synonyms	Tubulin alpha; Alpha Tubulin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1: 5000 - 1: 20000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of C.elegans cells using Tubulin alpha rabbit pAb(Caenorhabditis elegans). Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com