

## Tubulin beta Rabbit Polyclonal Antibody(Caenorhabditis elegans)

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP13266
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	20ul, 50ul, 100ul
<b>Immunogen</b>	Purified recombinant Full length Tubulin beta (Caenorhabditis elegans)
<b>Mol wt</b>	N/A
<b>Species reactivity</b>	C.elegans
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Tubulin beta
<b>Synonyms</b>	Tubulin beta; Beta Tubulin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

N/A


### Recommended Dilution

WB: 1: 5000 - 1: 20000

ELISA: 1: 20000

Not yet tested in other applications.

### Images



Western Blot analysis of C.elegans cells using Tubulin beta rabbit pAb(Caenorhabditis elegans). Secondary antibody was diluted at 1:20000

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)