

U2AFL Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13319
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human U2AFL
Mol wt	52690
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1
Synonyms	U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 ;CCCH type zinc finger; RNA-binding motif and serine/arginine rich protein 1;U2;RNU2 small nuclear RNA auxiliary

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com