

UBR4 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13348
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human UBR4
Mol wt	570130
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	E3 ubiquitin-protein ligase UBR4
Synonyms	E3 ubiquitin-protein ligase UBR4 ;EC 6.3.2.-;600 kDa retinoblastoma protein-associated factor;N-recognin-4;Retinoblastoma-associated factor of 600 kDa;RBAF600;p600;Zinc finger UBR1-type pr

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com