

ULK1(Phospho-Ser757) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13372
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized phosho peptide around human ULK1 (Ser757)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	ULK1
Synonyms	ULK1 ;Ser757 ; Serine/threonine-protein kinase ULK1; EC 2.7.11.1; Autophagy-related protein 1 homolog; ATG1; hATG1; Unc-51-like kinase 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western Blot analysis of various cell lysis(provided by our customer). Primary Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:10000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com