

YAP(Phospho-Ser397) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13520
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human YAP (Ser397)
Mol wt	N/A
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	YAP
Synonyms	YAP ;Ser397 ; Yorkie homolog; 65 kDa Yes-associated protein; YAP65

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development| growth| repair| and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com