

YB1(Phospho-Ser102) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13540
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized phospho peptide around human YB1 (Ser102)
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	YB1
Synonyms	YB1 ;Ser102 ; Tryptase alpha/beta-1; Tryptase-1; EC 3.4.21.59; Tryptase I; Tryptase alpha-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Tryptases comprise a family of trypsin-like serine proteases| the peptidase family S1. Tryptases are enzymatically active only as heparin-stabilized tetramers| and they are resistant to all known endogenous proteinase inhibitors. Several tryptase genes are clustered on chromosome 16p13.3. These genes are characterized by several distinct features. They have a highly conserved 3' UTR and contain tandem repeat sequences at the 5' flank and 3' UTR which are thought to play a role in regulation of the mRNA stability. These genes have an intron immediately upstream of the initiator Met codon| which separates the site of transcription initiation from protein coding sequence. This feature is characteristic of tryptases but is unusual in other genes. The alleles of this gene exhibit an unusual amount of sequence variation| such that the alleles were once thought to represent two separate gene

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year