

YY1 (Acetyl Lys351) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13583
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human YY1 (Acetyl Lys351)
Mol wt	45540
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	YY1
Synonyms	YY1 ;Acetyl Lys351; Transcriptional repressor protein YY1; Delta transcription factor; INO80 complex subunit S; NF-E1; Yin and yang 1; YY-1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

YY1 is a ubiquitously distributed transcription factor belonging to the GLI-Kruppel class of zinc finger proteins. The protein is involved in repressing and activating a diverse number of promoters. YY1 may direct histone deacetylases and histone acetyltransferases to a promoter in order to activate or repress the promoter| thus implicating histone modification in the function of YY1.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com