

ZDHC1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13655
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ZDHC1
Mol wt	53350
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Probable palmitoyltransferase ZDHHC1
Synonyms	Probable palmitoyltransferase ZDHHC1 ;EC 2.3.1.-;DHHC domain-containing cysteine-rich protein 1;Zinc finger DHHC domain-containing protein 1;DHHC-1;Zinc finger protein 377

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com